




Helping those who help patients

**ORDERING IS EASY!**  
 Online: [www.PathSUPPLY.com](http://www.PathSUPPLY.com)  
 Email: [Orders@PathSUPPLY.com](mailto:Orders@PathSUPPLY.com)



## Instrument Filters



Item #	Filter	Item #	Filter
<b>Filter-AF002</b> Replaces Part# FR-4	<b>AirFiltronix Ductless Fume Hoods</b> Activated Carbon 17"x10"x1 7/8" Gasket one side & particulate filter both sides	<b>FILTER-SA001</b> Replaces Part# 4349	<b>Sakura DRS-601 Stainer</b> Activated Carbon 25"x3 1/2"x3/8" gasket 1/8"
<b>FILTER-FG001</b> Replaces Part# 906	<b>Fume-Guard</b> Activated Carbon 11-1/4"x5-3/4"x1- 1/4"	<b>FILTER-SA002</b> Replaces Part# 7107	<b>Sakura EXPRESS X120</b> Activated Carbon 11-1/2"x9"x1 7/8"
<b>FILTER-LE001</b> Replaces Part# 045625623 & 14047432273	<b>Leica Autostainers ST4040, ST5010 &amp; ST5020</b> Activated Carbon 19"x4"x1-3/8" Gasket & prefilter	<b>FILTER-SA003</b> Replaces Part# 4899	<b>Sakura VIP</b> Activated Carbon 16 1/2"x8 1/2"x2 5/16" Cardboard Frame
<b>FILTER-LE002</b> Replaces Part# 046127121	<b>Leica TP1050</b> Activated Carbon 11 7/16"x5 5/8"x2"	<b>FILTER-SA004</b> Replaces Part# 4775	<b>SCA Coverslipper</b> Activated Carbon 8"x5"x1 3/8" gasket both sides
<b>FILTER-LE003</b> Replaces Part# 047634150 & 14047634150	<b>Leica ASP300</b> Activated Carbon 11 1/4"x5 3/4"x3/4"	<b>FILTER-SA005</b> Replaces Part# 2008	<b>Sakura DRS 2000 Stainer</b> Activated Carbon 24 1/2"x3 1/2"x1"
<b>FILTER-LE004</b> Replaces Part# 046427756	<b>Leica CV5000 Coverslipper</b> Activated Carbon 5 3/4"x5 3/4"x3/8"	<b>FILTER-SA006</b> Replaces Part# 6160	<b>Tissue Tek VIP6, Prisma Stainer, G2 Glass Cover Slipper, &amp; Express X50</b> Activated Carbon 14"x2"x4" 2 filters/package
<b>FILTER-LE005</b> Replaces Part# 042230673 & 14042230673	<b>Leica CV5030 Coverslipper</b> Activated Carbon 8"x7 1/2"x3/4"	<b>FILTER-SH001A</b> Thermo Shandon Item #990610	<b>Shandon Excelsior</b> CARBON FILTER – GLX Stainer, Hyperclean, Excel- sior, Varistain etc. Size 145 x 35 x 285 mm 5.71 " x 1.37" x 11.22"
<b>FILTER-LE006</b> Replaces Part# 0260434 & 526.0434	<b>Peloris II</b> Activated Carbon 15 3/4"x7 3/8"x1 1/4"	<b>FILTER-SH001B</b> Thermo Shandon Item #990612	<b>Shandon™ Excelsior™</b> PERMANGANATE FILTER - Hyperclean, Excelsior etc Size 145 x 35 x 285 mm 5.71" x 1.37" x 11.22"
<b>FILTER-LE007</b> Replaces Part# 14049543860	<b>Leica ASP6025 Automated Vacuum Tissue Processor 9.5" x 6.75" x 3/4" c/w gasket one side</b> Particulate filter both sides & pull strap	<b>FILTER-SP001</b> Replaces Part# 380043962	<b>Surgipath XYL-X50B Xylene Air Filter System</b> Activated Carbon 11 5/8"x3 7/8"x5 1/4" 3/4 hole on both ends with gasket
<b>FILTER-LE008</b> Replaces 042230674	<b>Formaldehyde Filter for Leica TP1020 Tissue Processor</b>		

Order online: [www.pathsupply.com](http://www.pathsupply.com) Order via email: [orders@pathsupply.com](mailto:orders@pathsupply.com) Order by phone: 800-631-3556